Urgent Medical Device Removal - Immediate Action Required

Imager™ 5F II Angiographic Catheter

February 11, 2020

Dear Vascular Surgery Department:

Boston Scientific Corporation (BSC) is initiating a removal of specific lots/batches of Imager II 5F Angiographic Catheters. BSC has noted an increase in the rate of tip detachment complaints involving units within these lots/batches. The preliminary investigation indicates that these batches meet design and manufacturing requirements, however external factors may have contributed to the tips of devices in these batches becoming brittle, leading to the tip detaching. No other Imager II Catheters are impacted by this removal.

The most common injury would be related to the tip detaching inside the patient resulting in either the need for intervention to retrieve the fragment or in the fragment remaining in the patient's vessel potentially requiring additional intervention and/or prolonged hospitalization. There is a possibility that a potentially life-threatening embolism of the device fragment could result.

This removal affects only the lots/batches listed in the table below. **Further distribution or use of any remaining product affected by this removal should cease immediately.**

Affected Product Listing

Product Description	Outer Package UPN #	Inner Package UPN #	GTIN	Lot/Batch #	Expiration Date
Imager™ II Angiographic Catheter	M001314051	M001314050	08714729354871	134092	23-Aug-2020
	M001314051	M001314050	08714729354871	134600	12-Sep-2020
	M001314061	M001314060	08714729354888	134011	20-Aug-2020
	M001314141	M001314140	08714729354963	133737	10-Aug-2020
	M001314341	M001314340	08714729355168	139512	12-Mar-2021
	M001314581	M001314580	08714729355403	134631	13-Sep-2020
	M001314591	M001314590	08714729355410	132447	13-Jun-2020
	M001314661	M001314660	08714729355489	132355	8-Jun-2020
	M001315151	M001315150	08714729355892	132823	26-Jun-2020
	M001315151	M001315150	08714729355892	133447	13-Jul-2020
	M001315151	M001315150	08714729355892	133448	16-Jul-2020
	M001315151	M001315150	08714729355892	134946	25-Sep-2020

Please ensure Imager II 5F Angiographic Catheters are stored per the DFU recommendations: Imager II 5F Angiographic Catheters must be stored in a cool, dry, dark place.

PLEASE NOTE: We are aware that hospitals often remove products from the outer carton and store them on shelves in the inner-pouch only. If this is a practice at your facility, it is **very important that you carefully use the product table and consider both the inner and outer packaging UPN codes when searching for affected product, as the UPN numbers on the inner and outer labeling are different. The product information listed on your specific Verification Form (enclosed with this letter) provides outer package product coding only** and should be utilized when reporting product to return.

Verify by product lot/batch number in the product table whether any product within your inventory is affected. If so, indicate on your Verification Form the quantity of units from each lot/batch that you will be returning. **As the product within these lots/batches are sold as 5-packs, it is important that all reported quantities represent the actual number of single units being returned and not the number of cartons/boxes or multi-packs.**

If you identify any product from the affected lots/batches within your inventory, please segregate the product immediately and return it to BSC in accordance with the enclosed instructions. **If you are a distributor, please note that the depth is to the hospital level and this notification should be forwarded to your customers. If you are a facility that has sent products to another hospital within your network, please ensure that this notification is forwarded to them.**

BSC is notifying regulatory authorities of this removal as required.

Please read carefully through the removal instructions included with this notification. Your local Sales Representative can answer any questions that you may have regarding this notification.

We regret any inconvenience that this action may cause, and we appreciate your understanding as we take action to ensure patient safety and customer satisfaction. We are committed to continuing to offer products that meet the highest quality standards that you expect from Boston Scientific Corporation.

Sincerely,

Nicole Maloney
Boston Scientific Quality Systems
763-494-1133
E-mail: BSCFieldActionCenter@bsci.com

Encl: Removal Instructions
Reply Verification Tracking Form
Field Action Return Shipping Label

Health care professionals and consumers may report serious adverse events (side effects) or product quality problems with the use of this product to Boston Scientific by calling 1-800-811-3211 and to the FDA's MedWatch Adverse Event Reporting program either online, by regular mail, fax or phone.

Online: www.fda.gov/MedWatch/report.htm

Regular Mail: use postage-paid FDA form 3500 available at www.fda.gov/MedWatch/getforms.htm and mail to MedWatch, 5600 Fishers Lane, Rockville, MD, 20852-9787

Fax: (800) FDA-0178

Phone: (800) FDA-1088

Removal Instructions

The Reply Verification Tracking Form enclosed with this letter must be completed and returned **even if you do not have any of the units from the recalled lot/batches.**

1. **Immediately discontinue use of and segregate recalled product.**
 - Immediately remove all affected recalled product from your inventory.
 - Segregate this product in a secure location for return to BSC.
2. **You must complete and email or fax your Reply Verification Tracking Form(s) immediately to get an RGA # before you return any product.**

- Indicate on your RVTF the quantity of units from each lot/batch that you will be returning.
- Return the RVTF as described below:

Email: BSCFieldActionCenter@bsci.com

or

Fax to: Field Action Center 1-866-213-1806

You will be contacted by BSC and provided a Returned Goods Authorization (RGA) Number within 24-48 hours **after your RVTF is received**

- **Report all inventory to be returned and report in single-unit quantity(ies), not boxes.**

Credits will be issued for all recalled product that is properly returned to BSC.

3. **Package/Ship the Recalled Product.**

- Package any product that is being returned in an appropriate shipping box.
- Affix the enclosed (red/white) shipping label to the outside of the shipping box.
- Write the **RGA number** in large print on the outside of the box, on the shipping label.
- Feel free to use our Federal Express Number 9205-2515-6 to return this package via second day delivery.
- Seal the box, and return to: **Boston Scientific Corporation**

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